

Bio Skin® Shinsplints Wrap

Application Instructions



1. Place sewn-in pad (end opposite the label) over painful area on the front of the shin.



2. While holding the pad end in place, wrap the support around the leg and secure the first hook strip onto the wrap.



3. Continue wrapping and tension the wrap to comfort and fasten end to secure.



Easy Care Instructions

Wash Bio Skin® regularly to prevent body salt and oil build-up that naturally occurs after the product has been worn a few times. Attach the hook to the loop to prevent lint build-up. This will give your Bio Skin® a longer life span.



Hand Wash



Hang Dry



Do Not Bleach



Do Not Dry Clean

Cropper Medical, Inc.
custserv@bioskin.com
800.541.2455
www.bioskin.com

